

# MATTINGLY LOW VISION TEST CARD

**Note:** The patient should use this card with their general wear multi-focal spectacles (usually bifocals). The patient should progress down the card to the smallest line that is comfortable to read. The approximate magnification needed to read small print (1M) is given to the right of the line read. The magnifier icon indicates possible optical solutions such as magnifiers and spectacles, and the video reader icon indicates a possible need for this type of device.

Reading



(10M) 10X/40D

Requires



(8M) 8X/32D

Adequate



(6M) 6X/24D

Lighting and



(5M) 5X/20D

Magnification.



(4M) 4X/16D

A magnifier should



(3M) 3X/12D

Allow you to read your mail.



(2M) 2X/8D

This print is similar to small headlines.



(1.5M) 1.5X/6D

Congratulations! You are now reading small print.

(1.M)